MANUAL DEL USUARIO
OWNER'S MANUAL

KMINATM

MULETA KMINA PRO
KMINA PRO CRUTCH

CE

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www.kmina.es

+34 960 046 322

HOW TO ADJUST YOUR KMINA PRO CRUTCHES**1 ADJUST THE HEIGHT**

Straight shoulders.

Shoulders are too far down.

Shoulders are too far up.

2 ADJUST THE HANDLE POSITION**1°** Place the elbow at the suggested zone.**2°** Adjust it to the right position (see examples).**3°** Secure the position with the screw.

The handle position is correct.

The handle position is too close. The handle position is to far away.

***MAKE SURE THAT THE FASTENING STRAP IS PROPERLY PLACED**
We recommend to place the opening part of the fastening strap on the outside.

Right crutch.

Left crutch.

HOW TO USE YOUR KMINA PRO CRUTCHES

1 PLACE YOUR ELBOW CLOSE TO YOUR BODY

 YES NO

2 WALK WITH THE TIP 30 CM FROM YOUR BODY

 YES30 cm NO NO

3 WALK WITH YOUR ELBOWS CLOSE TO YOUR SIDES, RAISE YOUR HANDS UP AND MOVE THE TIP FORWARD

 YES NO

4 TAKE SHORT STEPS

 YES NO

5 TAKE OFF YOUR CRUTCH CORRECTLY

The flexible fastening strap allows you to take off your crutch from behind or from the side. For a better grip and greater safety you may experience some pressure.

HOW TO ADJUST. MAKE SURE THAT:

- Your elbow is placed in the suggested zone.
- The handle padding is well adjusted and does not slip out of place.
- Nothing is unadjusted or uneven and there are no cracks in the metal.
- All springs and screws are tight and well regulated.
- The forearm support is at a 30 degree angle with respect to the horizontal axis.
- The grip of the handle is well regulated, placing the elbow in the indicated area.
- The maximum forearm length for the crutches should not exceed 24cm. The minimum forearm length requirement is 20cm. Measure from the elbow to the wrist.

HOW TO USE. MAKE SURE THAT:

- You do not have any disability or are under the influence of any medicine, etc., that a medical professional or therapist would deem to be in conflict with the safe use of the product.
- You have consulted a doctor or therapist and have been informed of the proper and safe way to walk with these crutches.
- You have correctly adjusted the height of the crutch, which has 17 adjustments every 1.5cm. The minimum height of the crutch is 96cm and the maximum height is 121 cm.
- If you use two crutches, ensure they are both set to the same height and all springs are working correctly.
- Body weight should be evenly distributed throughout your forearm and not on the hand.
- In the event of a sudden accident or shock, do not hold the crutches.
- Release your hands and try to place them on the ground.
- The maximum weight capacity for this crutch should not exceed 100kg.
- The maximum height capacity for this crutch should not exceed 190cm.
- The minimum height requirement is 150cm.

HOW TO USE. MAKE SURE THAT:

- Try and avoid potential risks and hazards that would put your safety in jeopardy eg, slippery surfaces, irregular surfaces and objects on the ground.
- Using these crutches does not eliminate the potential risk of falling.
- Whilst these crutches have been designed to be more comfortable for the patient, careful use should be carried out at all times.
- These crutches should be kept clean at all times.

**IMPROVING MOBILITY DEVELOPMENT SL.**

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