



SYST'AM® P951D / BACK SUPPORT & THERAPEUTIC COMBI-PACKS

BACK SUPPORT MADE OF VISCOELASTIC FOAM WITH MEMORY EFFECT AND MULTI-BEARING SURFACE

MATERIALS

80
kg/m³



Visco foam



Removable
POLYMAILLE®
cover

Foam maintenance:



(Do not immerse in water)

Cleaning of the cover:



INDICATIONS

P951D / Back support and Therapeutic Combi-Packs can be used to:

- 1- Improve positioning in wheelchair.
- 2- Improve seating orientation
- 3- Contribute reducing forward sliding.
- 4- Contribute releasing back pains.



UNIVERSAL ATTACHMENT SYSTEM

→ Enables adaptation of the back support to most wheelchair models

AVAILABLE VERSIONS



P951D / Back Support
With POLYMAILLE® Cover.



Therapeutic Combi-Packs
(See P. 46)

DESIGNATION	SIZE	ITEM CODE	SIZES W x L x H (cm x inches)
P951D / Back support Small with POLYMAILLE® Cover	Small	P951DS1HW	44 x 39 x 9 cm / 17 x 15 x 3,5"
	Standard	P951DIHW	49 x 44 x 9 cm / 19 x 17 x 3,5"
	Large	P951DL1HW	49 x 50 x 10 cm / 19 x 20 x 4"

